



Simply The ONE



Reliable

Universal

Aesthetic

Technique insensitive*

Simple

G-CEM ONE

Self-adhesive resin cement

* Miyazaki S (2020). Excellent curing properties: chemical polymerization and comfortable clinical practice compatibility of G-Cem ONE. Clinical Advance; 7:167-169.



Since 1921
100 years of Quality in Dental

The ONE... that simplifies!

G-CEM ONE is a truly **universal, non-technique sensitive** cement with **high bond strength** and excellent self-curing ability. It simplifies **all your cementation procedures** and gives you more **predictable results** and peace of mind. G-CEM ONE can be used for a **wide range of indications**, including challenging clinical situations thanks to the **GC Touch Cure**.

One GC Touch Cure

The chemical initiator in **Adhesive Enhancing Primer** will accelerate the chemical cure of the cement from the tooth surface to ensure stable bonding.




Reliable

Optimal bond strength for retentive and non-retentive¹ preparations thanks to the unique **GC Touch Cure** ability of the **optional Adhesive Enhancing Primer**.

Technique insensitive

Exceptional self curing[^] for your peace of mind, tack cure feature for easy excess removal and clean up, moisture and saliva tolerance¹, optimised flow of the paste and virtually no post-operative sensitivity³.

One optimised self-curing ability[^]



Source: Miyazaki S (2020). Excellent curing properties: chemical polymerization and comfortable clinical practice compatibility of G-Cem ONE. *Clinical Advance*; 7:167-169.

¹ With Adhesive Enhancing Primer. GC R&D Data on file.

² With G-Multi Primer. GC R&D Data on file.

³ Field survey results. Using G-CEM ONE only. GC Data on file.



The ONE... that simplifies!

Adhesive Enhancing Primer (Optional)

For all abutment tooth and cavities

- MDP: bond to tooth structure, zirconia and non-precious metal
- 4-MET: bond to tooth structure
- MDTP: bond to precious metal
- Accelerator → Strong bonding to abutment tooth and cavity by **GC Touch Cure** feature.

No light curing, simply drying.

Wait 10 sec

Strong air 5 sec



Universal

Excellent bond strength to tooth and to all substrates:

Enamel, Dentin, Zirconia, Metal, Glass ceramics², Hybrid ceramics².

Aesthetic

Invisible, wear resistant margin for aesthetic outcomes.

Simple

Reduced inventory for reproducibility and optimised workflow.

One performance
regardless of the translucency



Metal

Images courtesy of
Dr. Kazunori Otani, Japan





Lithium disilicate


Images courtesy of
Dr. Akimasa Tsujimoto, Japan

The ONE... you can trust!


One cement for superior bond strength to all substrates



One primer for clinical challenges




Immediate bond strength to avoid early debonding in non-retentive cases and CAD/CAM milled crowns, for the comfort of your patients.



The Adhesive Enhancing Primer (AEP) ensures high bond strength, even in presence of saliva contamination.

One invisible margin for aesthetic outcome



Low film thickness (5 µm) and high wear resistance keep the invisible cement line.



Before treatment



After six months

Images courtesy of Dr. Masayuki Okawa, Japan



With optional Adhesive Enhancing Primer

The ONE... clinical partner!

One optimised flow



One easy excess removal with 1 second tack curing



Images courtesy of
Dr. Kazunori Otani, Japan

One for all Clinical case scenarios

Posterior Crown (Zirconia)



Prepared tooth cleaned and dried.



Restoration cemented.



Final result.

Images courtesy of
Dr. Kazunori Otani, Japan



Anterior Crown (Lithium Disilicate)



Prepared tooth cleaned and dried.



Application of Adhesive Enhancing Primer.




Final result.

Images courtesy of
Dr. Kazunori Otani, Japan




Posterior Onlay (Lithium Disilicate)



Initial situation.



Restoration cemented with use of Adhesive Enhancing Primer.



Final result.

Images courtesy of
Dr. Riccardo Tempesta, Italy



With optional Adhesive Enhancing Primer

Packaging

Description

G-CEM ONE Starter Kit A2

1 x G-CEM ONE A2 4.6g (2.7mL)
1 x G-CEM ONE ADHESIVE ENHANCING PRIMER 2mL
8 x G-CEM Automix Tip Regular
2 x G-CEM Automix Tip Endo

G-CEM ONE Starter Kit Translucent

1 x G-CEM ONE Translucent 4.6g (2.7mL)
1 x G-CEM ONE ADHESIVE ENHANCING PRIMER 2mL
8 x G-CEM Automix Tip Regular
2 x G-CEM Automix Tip Endo

G-CEM ONE Twin Refill A2

2 x G-CEM ONE A2 4.6g (2.7mL)
15 x G-CEM Automix Tip Regular
5 x G-CEM Automix Tip Endo

G-CEM ONE Twin Refill Translucent

2 x G-CEM ONE Translucent 4.6g (2.7mL)
15 x G-CEM Automix Tip Regular
5 x G-CEM Automix Tip Endo

G-CEM ONE Single Refill AO3

1 x G-CEM ONE AO3 4.6g (2.7mL)
8 x G-CEM Automix Tip Regular
2 x G-CEM Automix Tip Endo

G-CEM ONE Single Refill White Opaque

1 x G-CEM ONE White Opaque 4.6g (2.7mL)
8 x G-CEM Automix Tip Regular
2 x G-CEM Automix Tip Endo

G-CEM ONE ADHESIVE ENHANCING PRIMER, 4mL



G-CEM ONE Starter Kit




G-CEM ONE Refill



G-CEM ONE ADHESIVE ENHANCING PRIMER

Shade range



Related products



G-Multi PRIMER
Universal Primer



FREEGENOL
Temporary Luting Cement



Fuji TEMP LT
Long Term Temporary Cement

G Multi PRIMER, FREEGENOL, Fuji TEMP LT are trademarks of GC.



Since 1921
Towards Century of Health

GC Australasia Dental Pty Ltd
1753 Botany Road Banksmeadow
NSW 2019 Australia

T: +61 2 9301 8200
E: info.australasia@gc.dental
www.gcaustralasia.com