



GAUZE SWAB - 8PLY




100% cotton and are French folded to reduce fraying and linting. The gauze swabs are sterile and 8ply in thickness and ideal to be used as a wound dressing or absorbing fluid.

PRODUCT FEATURES

- o Highly absorbent

SUITABLE FOR

- o Cleaning wounds
- o Sensitive skin
- o Soak up fluid
- o Padding
- o Wound protection barrier



REORDER	PCS	SIZE	COLOUR	PCS/UNIT	UNITS/CTN
GSG2425	5pcs	10cm x 10cm	Green	50	9
GSR2223	3pcs	5cm x 5cm	White	50	30
GSR2225	5pcs	5cm x 5cm	White	50	24
GSR2322	2pcs	7.5cm x 7.5cm	White	50	24
GSR2323	3pcs	7.5cm x 7.5cm	White	50	24
GSR2325	5pcs	7.5cm x 7.5cm	White	50	24
GSR2422	2pcs	10cm x 10cm	White	50	16
GSR2423	3pcs	10cm x 10cm	White	50	8
GSR2425	5pcs	10cm x 10cm	White	50	8

All product information is correct at date of publication

