

SAFETY DATA SHEET

Section 1: IDENTIFICATION

INSTRUMAX PINK

Recommended Use: Instrument Disinfectant - Low Level.

Product Code: 090041 (2x5L), 090043 (6x500mL Trigger Spray), 090045 (12x500mL Squeeze).

Whiteley Medical

A division of Whiteley Corporation Pty Ltd (A.C.N. 000 906 678)
Postal Address: P. O. Box 1076 North Sydney NSW 2059
Telephone Number: (02) 9929 9155 Facsimile: (02) 9929 9077

Web: www.whiteley.com.au

Emergency Telephone Number: Poisons Information Centre (National) 131126

Distributed in New Zealand by: Whiteley Corporation NZ Limited

Address: PO Box 22-519 Otahuhu, Auckland, New Zealand 1640
Telephone Number: 092737313, 0800 257 352
Emergency Telephone Number: New Zealand Poisons Centre 0800 POISON or 0800 764 766

Section 2: HAZARDS

GHS Classification

Not a dangerous substance or mixture according to GHS or HSNO criteria.

Signal Word None

Hazard Statements

None

Precautionary Statements

None

Response Statements

None

Storage Statements

None

Disposal Statements

None

HSNO Classifications

N/A

Issue No. 2 GHS Page 1 of 5

SAFETY DATA SHEET INSTRUMAX PINK

Section 3: COMPOSITION INFORMATION

Ingredient CAS No **Proportion** Ethanol 64-17-5 <10% Ingredients deemed not to be hazardous Not applicable To 100%

Section 4: FIRST AID

Eye (Contact) Hold eyelids apart and flush the eye continuously with running water.

Skin (Contact) Remove contaminated clothing and flush skin and hair with running water.

Inhalation(Breathing) Remove to fresh air.

Ingestion (Swallowing) DO NOT induce vomiting. Give water to drink if conscious. For advice,

contact a Poisons Information Centre (Phone 131126) or a doctor.

Advice to Doctor Treat symptomatically for neutral detergent.

First Aid Facilities Ensure an eye wash and safety shower is available and ready for use.

Additional Information No aggravated medical conditions are known to be caused by exposure to

this product.

Section 5: FIREFIGHTING MEASURE

Suitable Extinguishing

Media

Use water fog (or if unavailable fine water spray), dry chemical, carbon

dioxide or alcohol stable foam.

Hazards From

Combustion Products

Carbon dioxide and carbon monoxide may be produced in the case of fire

or during thermal decomposition.

Precautions For Fire Fighters and Special Protective Equipment

Keep containers cool by spraying with water.

Additional Information Hazchem Code – Not applicable.

Section 6: ACCIDENTAL RELEASE MEASURES

Emergency Procedure SAA/SNZ HB76: Dangerous Goods - Initial Emergency Response

Guide - Not applicable.

Spills / Clean up Clean up personnel should wear full protective clothing. Restrict access

> until completion of clean up. Then ensure adequate ventilation. Stop leak if safe to do so. Contain spill with absorbent material, such as towelling, sand, vermiculite or other inert material. Prevent spill entering Collect and dispose of clean up stormwater drains or waterways. material according to local regulations. Wash away remnants with copious amounts of cold water to sewer. Clean area by working from the periphery to the centre of spill or from the edge of the room to the

centre.

Section 7: HANDLING AND STORAGE

Precautions for Safe

Handling

Contact Whiteley Corporation sales representative for advice when using this product for any application other than that outlined on the label

or technical bulletin.

Issue No. 2 GHS Page 2 of 5 Any non-intended or non-authorised use of this product may result in personal injury or damage to equipment.

personal injury or damage to equipment.

Store product in original container.

Wash thoroughly after handling product.

Conditions for Safe Storage

Store in a cool, dry, well ventilated area. Keep container tightly sealed.

Store below 30°C out of direct sunlight.

Section 8: EXPOSURE CONTROL/PERSONAL PROTECTION

National Exposure Standards – Source: National Exposure Standards for Atmospheric Contaminants in the Occupational Environment [NOHSC:1003] & New Zealand Workplace Exposure Standards 2002.

IngredientCAS NoES-TWAES-STELEthanol64-17-51000ppm-

Biological Limit Values No data available.

Engineering Controls Use only in a well ventilated area.

Personal Protective
Equipment

Eye/face protection – Safety glasses / face shield / chemical resistant goggles should be worn to prevent eye contact.

Skin protection – Use Nitrile gloves or similar to prevent skin contact.

Respiratory protection – Respirator is not usually necessary but if required use a half face filter respirator suitable for organic vapours.

Section 9: PHYSICAL AND CHEMICAL PROPERTIES			
Appearance	Clear pink liquid	Odour	Citrus
Odour Threshold	No data available	рН	7.00 - 7.40
Freezing Point	Approx. 0°C	Initial Boiling Point	Approx. 100°C
Boiling Range	No data available	Flash point	No data available
Evaporation Rate	No data available	Flammability	No data available
Upper flammability limit	No data available	Lower flammability limit	No data available
Vapour Pressure	No data available	Vapour Density	No data available
Relative Density	0.985 – 0.995	Solubility	Completely miscible with water
Partition Coefficient: n-octanol/water	No data available	Autoignition temperature	No data available
Decomposition temperature	No data available	Viscosity	Approx. 1 cPs

	Section 10: STABILITY AND REACTIVITY	
Chemical Stability	Product is stable and will not undergo any hazardous reactions under normal conditions of use and storage.	
Conditions to avoid	Avoid high temperatures (store below 30°C) Do NOT store in direct sunlight.	

Issue No. 2 GHS Page 3 of 5

SAFETY DATA SHEET

INSTRUMAX PINK

Incompatible materials None known. Do not mix with other chemicals.

Hazardous decomposition products

None known.

Hazardous reactions None known.

Section 11: TOXICOLOGICAL INFORMATION

No data available. **Acute toxicity**

Skin corrosion/

irritation

No data available.

Serious eye damage/

irritation

May cause irritation.

Respiratory or skin

sensitisation

May cause irritation.

Germ Cell Mutagenicity

No data available.

Carcinogenicity

No components of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen

by IARC.

Reproductive toxicity No data available.

Specific Target Organ

Toxicity - Single **Exposure**

No data available.

Specific Target Organ Toxicity - Repeated

Exposure

No data available.

Aspiration Hazard

No data available.

Section 12: ECOLOGICAL INFORMATION

Ecotoxicity No data available.

Persistence and degradability

No data available.

Bioaccumulative

potential

No data available.

Mobility in soil No data available.

Other adverse effects No data available.

Section 13: DISPOSAL CONSIDERATIONS

Disposal method Disposal to sewer is normally recommended with copious amounts of

water.

Refer to State/Territory Land Waste Management Authorities if applicable. Containers are recyclable and can be disposed of by a licensed waste contractor. Containers can be disposed of to general

Issue No. 2 GHS Page 4 of 5 waste or rinsed thoroughly and recycled.

Special precautions Suitable for incineration by approved agent.

Section 14: TRANSPORT INFORMATION

Not classified as Dangerous Goods by the criteria of the Australian Dangerous Goods Code (ADG

Code).

UN Number Not applicable. **UN Proper Shipping** Not applicable.

Name

Class and subsidiary

Not applicable.

risk

Packing Group Not applicable.

Special precautions for Not applicable.

user

Hazchem Code Not applicable.

Section 15: REGULATORY INFORMATION

Poisons Schedule (SUSDP): Not applicable. HSNO Approval Code: Not applicable.

All ingredients are listed in the Australia Inventory of Chemical Substances (AICS).

This document has been produced in accordance with the requirements of the Globally Harmonised System of Classification and Labelling.

Section 16: OTHER INFORMATION

Date of

preparation: 12th June 2018

© This document remains property of Whiteley Corporation Pty Ltd.

Alterations are not permitted without prior written authorisation from Whiteley

Corporation Pty Ltd.

Issue No. 2 GHS Page 5 of 5