


Your practice is *our* inspiration.™

**Take 1 Advanced** | Designed for the way you work



# Take 1 Advanced

A perfect fit for you, and your patients




Our expertise in high-performance impression materials, along with recent technological advancements, delivers a VPS material with the optimum combination of physical properties: strength, elasticity, dimensional stability, and the ability to register detail in any environment. Take 1 Advanced has the most extensive selection of set times, delivery systems and consistencies on the market today.

## Research

As with any dental material, the physical and mechanical properties of impression material dictate how well it will perform in different clinical situations. While some materials focus on a single attribute, it's critical for these properties to work well together. Take 1<sup>®</sup> Advanced<sup>™</sup> features individual properties that rate as good or better than any other VPS on the market today, with collective performance exceeding that of competitive products.

### Contact Angle

The patented formulation consisting of a polymerizable hydrophilic modifier and surfactants creates a highly wettable surface that breaks surface tension and easily displaces saliva and other fluids to capture outstanding detail in any environment.



Internal testing. Data available upon request.

## Features

**Outstanding tear strength and elongation.** Strong, flexible material stretches around undercuts without tearing or distorting.

**Hydrophilic wash material.** Displaces intraoral fluids to produce impressions with great detail and no voids or gaps.

**Widest selection of viscosities, set times and delivery methods.** Helps overcome common difficulties in impression taking, and matches doctor preference.




Internal testing. Data available upon request.

### Tear Strength

Take 1 Advanced Wash material resists tearing around undercuts upon removal from the mouth. That combined with a high elongation percentage ensures preparation margins remain intact, even in tough cases below the gumline.

## Material Set Times



### Options: Delivery, Viscosity and Set Times

LB Wash	RB Wash	Med/Mono	HB Tray	Rigid Tray	Putty	Bite Reg
Flowable, highly flexible and hydrophilic	Highly thixotropic with incredible tear strength	The perfect balance between flow and body	Outstanding dimensional stability and elasticity	Designed to provide stability & support closed bite trays	Ultra-high viscosity for extra mechanical compression	Exceptional resistance to deformation to protect occlusal integrity

	LB Wash	RB Wash	Med/Mono	HB Tray	Rigid Tray	Putty	Bite Reg
<b>Unidose®</b> Kerr's unique Unidose delivery offers dual benefits: targeted syringe delivery and the elimination of cross-contamination concerns	Super Fast	Super Fast					
	Fast Set	Fast Set					
	Regular Set	Regular Set					
<b>Cartridge</b> Convenient 50 ml cartridges make mixing & application of Take 1 Advanced precise and easy	Super Fast	Super Fast	Super Fast	Super Fast	Super Fast		Fast Set
	Fast Set	Fast Set	Fast Set	Fast Set	Fast Set		
	Regular Set	Regular Set	Regular Set	Regular Set			
<b>Volume</b> For busy offices, Take 1 Advanced Volume delivery speeds up impression taking, reduces hand fatigue and lowers cost per use			Super Fast	Super Fast	Super Fast		
			Fast Set	Fast Set	Fast Set		
			Regular Set	Regular Set			
<b>Hand Mix</b> For instances when a higher viscosity tray material is desired and other mixing methods are unavailable						Regular Set	


Your practice is our inspiration.™



# Take 1 Advanced Bite Registration

Take 1 Advanced Bite completes the Take 1 Advanced portfolio providing you a complete solution within a single family of products.

Your practice is *our* inspiration.™



## Ordering

### 34064 Cartridge Intro Kit

Contains: 1 Cartridge (50 ml) LB Wash (Fast Set)  
 1 Cartridge (50 ml) RB Wash (Fast Set)  
 3 Cartridges (50 ml) Tray (Fast Set)  
 1 Cartridge (50 ml) Bite Registration  
 1 Extruder  
 6 Small Mix Tips  
 6 Intraoral Syringe Tips (Small Universal)  
 6 Large Mix Tips

### Refills

#### Two Packs: 2 Cartridges (50 ml), 6 Mixing Tips

33953 Light Body Wash FS  
 33951 Light Body Wash RS  
 33955 Light Body Wash SF  
 33971 Regular Body Wash FS  
 33969 Regular Body Wash RS  
 33973 Regular Body Wash SF  
 (Wash viscosity 2-packs also contain intraoral tips)

33959 Medium FS  
 33957 Medium RS  
 33961 Medium SF  
 33965 Tray FS  
 33963 Tray RS  
 33967 Tray SF  
 34457 Rigid Tray FS  
 34458 Rigid Tray SF  
 33975 Bite Registration

#### Value Packs: 24 Cartridges (50 ml)

33954 Light Body Wash FS  
 33952 Light Body Wash RS  
 33956 Light Body Wash SF  
 33972 Regular Body Wash FS  
 33970 Regular Body Wash RS  
 33974 Regular Body Wash SF  
 33960 Medium FS  
 33958 Medium RS  
 33962 Medium SF  
 33966 Tray FS  
 33964 Tray RS  
 33968 Tray SF  
 35186 Rigid Tray FS  
 35187 Rigid Tray SF

#### Unidose: 20 Syringes (2 ml)

33995 Light Body Wash FS 20 pack  
 33994 Light Body Wash RS 20 pack  
 33996 Light Body Wash SF 20 pack  
 33998 Regular Body Wash FS 20 pack  
 33997 Regular Body Wash RS 20 pack  
 33999 Regular Body Wash SF 20 pack

#### Volume: 2 Hard Cartridges (380 ml)

34695 Tray RS  
 34697 Tray SF  
 34696 Tray FS  
 34690 Medium/Monophase FS  
 34688 Medium/Monophase RS  
 34689 Medium/Monophase FS  
 34825 Rigid Tray FS  
 34826 Rigid Tray SF

#### Putty: 1 Base (400 g), 1 Catalyst (400 g)

34070 Putty

## Accessories

28747 Extruder  
 29363 Mix Tips, Small (pkg of 48)  
 29333 Mix Tips, Large (pkg of 48)  
 27517 Intraoral Syringe Tips, Small Universal (100 disposable)  
 20731 Intraoral Syringe Tips, Large (100 disposable)  
 16622 VPS Adhesive (15 ml bottle)  
 25777 VPS Adhesive (59 ml bottle)  
 35115 Volume Machine by Renfert  
 34746 Dynamic Mixing Tips for Volume (50 tips, 2 bayonet rings)

©2012 Kerr Corporation Lit. No. 34117D  
 All trademarks are property of their respective owners.