

PRIME LINE - PREMIUM RANGE

Your Silent Partner

• Only 55dB Operating Noise •



Power-Torque head (25W) and Small head (18W)



Premium spindle with carbide insert for high concentricity



Pressure reduction system in KaVo* multiplex fitting



Integrated non-return valve to prevent cross-infection



Ceramic bearings for longer life and quiet operation



40,000 LUX solid fibre optic cellular glass rod



Scratch resistance titanium coated prime-grip



Efficient cooling with 4-port or 3-port spray



~~\$1,265.00~~

25% OFF

Now Only \$889.45 Inc GST

while stock lasts

German handpieces at affordable prices and backed by factory warranty

LED
Compatible Handpieces

Available fittings:

KaVo* Coupling:
NSK* Coupling:
W&H* Coupling:
Sirona* Coupling:

Power Head
25W

HP21KL
HP21NL
HP21WL
HP21SL

Small Head
18W

HP22KL
HP22NL
HP22WL
- NA -

Authorised Dealer:

No other offers can be used with this deal.



MATRIX
DENTAL SERVICES

• Specials available while stocks last •

www.mk-dent.com.au

* KaVo, Sirona, NSK, and W&H are registered trademarks of KaVo Dental GmbH, Sirona Dental Systems GmbH, NSK, and W&H Dental Work Bürmoos Inc. respectively. MK-dent is in no way affiliated with these companies.

ULTRASONIC SCALER

• Perfect Cleaning •

Get 1 DLC coated black tip FREE (similar to EMS[®] A-type) with each Scaler



For EMS[®] connection

~ SC21E (SETSC2) \$330.00



EMS[®] connection
Non-Optic, type EN-041

~ SC21E2 (SETSC5) \$402.00



EMS[®] No Pain connection
Non-Optic, type EN-061

~ SC21EL (SETSC1) \$505.00



EMS[®] No Pain connection
with LED, type EN-060A

23 different types of tips available
for EMS[®] & Satelec[®] Scalers
Please visit www.mk-dent.com

For Satelec[®] connection

~ SC21S (SETSC4) \$330.00



Acteon-Satelec[®] connection
Non-Optic, type F 12 281

~ SC21SL (SETSC3) \$505.00



Acteon-Satelec[®] connection
with LED, type F 12 609

For KaVo[®] connection

~ SC21KL* (SETSC6) \$505.00



KaVo[®] connection
with LED, Ref 1.007.3995

*Only for use with MK-dent STxxE/ED
or EMS[®] original scaler tips

Get 3 Proply Powder FREE with each Proply Unit



Now Available at a
REDUCED PRICE

\$1,082.00
Inc GST

- ▶ Ergonomically shaped, Compact and Light-Weight Air-Driven Polishing Device.
- ▶ **360°** rotatable outer handle due to Unique **RotaGrip[®]** for Simple and Fatigue-Proof Operation.

Available coupling fittings:

KaVo* (PR1011K), W&H* (PR1011W),
Sirona* (PR1011S), NSK* (PR1011N),
BienAir* (PR1011B)



PP1011
Mint



PP1012
Lemon



PP1013
Strawberry

Extra Fine Grain (12-14 μm)

\$60.50 each Inc GST

ECO LINE

• High Flexibility - Low Price •



40,000 LUX solid fibre optic cellular glass rod



Power head (21W) and Small head (16W) options



Ceramic bearings for longer life and quiet operation



Efficient cooling with 4-port or 3-port spray



Smooth grip with modern and ergonomic design

**BUY 3+
SAVE
\$195.00**

**BUY 3+
\$543.00
each**

**BUY 1+
\$608.00
each**



<u>Available for</u>	<u>Power Head</u>	<u>Small Head</u>
KaVo* Coupling:	HE21KL	HE22KL
NSK* Coupling:	HE21NL	HE22NL
W&H* Coupling:	HE21WL	HE22WL
Sirona* Coupling:	HE21SL	- NA -

• Specials available while stocks last •

Buy 4 - Get 1 FREE (Mix & Match)

BASIC LINE

• Economical Option •



LB11/11

- 20,000 rpm
- Blue (1:1)
- Push Button

\$360.00



LB41/11P

- 20,000 rpm
- Green (4:1)
- For Snap-On Prophylaxis Cups

\$360.00



LB41/11

- 20,000 rpm
- Green (4:1)
- Push Button

\$412.00



LB01

- 20,000 rpm
- Blue (1:1)

\$360.00



MONT BLANC®

10000XL - Red Band (1:5)
Contra angle with Fibre Optic

Special Price

\$1,275.00 Inc GST

While Stock lasts

RRP: \$1,195.00 Inc. GST

Save 20%



Quadruple (4 hole)
Water Spray outlets



Anti Particle Filter
(Patent Pending)



25,000 L UX Cellular
Glass Rod Optics



Small Head with
Ceramic Bearings



i.coat® or
Ideal Handgrip




Lightweight due to
Titanium components




Fits ALL Electric &
Air (e-type) Motors




3 Easy Steps



1. Please select the
recommended Torque



2. Start rotating



3. Screw automatically
stops when selected
torque is reached

Lightweight
ONLY 135 gms




Made in France

Model:15501
\$917.00 Inc GST

Universal Torque Wrench

TORQ CONTROL®

Precision instrument with a Torque Adjustment System for tightening of prosthetic components at a specific torque, for any implant system



Automatic declutching when
selected torque is reached

- Compliance with Protocols for screw tightening
- Limits the risk of screw fracturing or loosening




Universal connection of
prosthetic mandrels

- 7 different Torque levels (from 10 to 35 N.cm)
- Can be used with ALL Implant Systems



Micro-Head and Ergonomic
handle

- Easy access in restricted area with Excellent Visibility



Easy and Quick to use

- Reduced motion compared with using a ratchet wrench while tightening