

pspix2

By your side, to care for everyone

ACTEON® place its 20 years of expertise at your service and those of your patients.

The first simple and intuitive personal imaging plate scanner provides outstanding images and was designed to free practitioners from constraints.

Whatever your patient's morphology or pathology, the PSPIX² adapts to your needs. As your personal assistant or in collaborative mode, experience high productivity without compromising on quality.

- 1 - **Striking Contrast** for a more reliable Diagnosis
- 2 - **The Smallest Scanner** on the market
- 3 - **The Fastest Scanner** on the market
- 4 - The **most intuitive** Scanner
- 5 - **Thin and flexible imaging plates** of various sizes for adaptability and outstanding patient comfort.



Special Offer Buy PSPIX2 \$12,990.00 RRP

x-mind unity

Versatility for "extended" possibilities

With its elegant style coupled to the cutting edge patented ACE® technology (Automatic Exposure Control) and the possibility to integrate SOPIX® Plug-In sensor, X-MIND® Unity brings the standards of protection, accuracy and comfort to a next level.

X-Mind® unity is available in different configurations allowing to fulfill any operational need. Arm lengths: 0.40 m, 0.80 m, 1.10 m - installation type: top wall, bottom wall, mobile stand.

\$7,276.00 RRP



BUY 1 x Plug N Sensors (Mix N Match Ok) GET 1 x Sopix 2 Plug In Dock FREE

x-mind ac x-mind dc

acteon

Experience intuitive excellence

Designed to simplify daily operations, prioritise patient safety and enhanced dental workflows X-MIND AC/DC is the perfect choice for modern dental professionals.

Three arms lengths are available : 41 cm (16.2"), 82.5 cm (32.5") and 110 cm (43.5") and can be: Top wall mounted Bottom wall mounted Dental unit, ceiling and mobile version are also available.

X-MIND DC available with **Remote Wireless Control**

- No more wiring issues
- Cost effective
- Enhanced Convenience



X-Mind DC \$5,199.00 RRP New X-Mind DC Wireless Button \$945.00 RRP

C50 Full HD Camera
with diagnostic Aid Option

The **C50 Camera** with diagnostic-aid option is the only camera on the market to cover so many pathologies, featuring 4 integrated vision modes for seamless navigation.

Remarkable clinical vision

- **Full HD resolution** for superior image quality
- **Autofocus** for a clear image
- **4 camera shots available;** macrovision, single tooth, smile and portait

Remarakable clinical versatility



<p>DAYLIGHT MODE</p>	Record of natural colors for a better interpretation of tooth shade, pathology identification and patient follow-up	<p>DAYLIGHT+ MODE</p>	High contrasts emphasizing anatomic fissures and details of the dentin structure/depth.
<p>PERIO MODE</p>	Enhancement of dental plaque and gingival inflammation thanks to chromatic amplification.	<p>CARIO MODE</p>	Caries detection by fluorescence technique by LEDs.

To order: C50 with diagnostic Aid S_C50_0001 C50 without diagnostic Aid S_C50_0002

For more Information or to place orders, call **Acteon ANZ**
Toll Free **1300 ACTEON** (1300 228 366) or Sydney (02) 9669 2292 or Email acteon.anz.promo@acteongroup.com

x-mind prime 2

acteon 
imaging suite

56 x 2D and 3D programs



- Easy interface to use
- Complete 3D diagnostic tools, powered by AIS3DAPP, to fully exploit the clinical content from your cases
- The Implant Planning Wizard combined with a Cloud implant library allows you to easily and safely plan your surgery in just a few clicks
- The Surgical Guide Design option brings your planning to the surgical room with the highest fidelity possible with a custom designed surgical stent
- Effective patient engagement: export to portable viewer and involve the patient in the diagnostic process and raise his awareness to boost case acceptance

Strong hardware combined with smart software to fast diagnostic.

Elevate your clinical performances

Whatever your specialty in dentistry, the quality of the diagnosis is a pillar of a clinical workflow to guarantee successful analysis, treatment plan, procedure and outcomes.

On X-MIND® prime, a complete choice of exam modalities combined with a premium-class image quality will get you the ability to successfully handle any patient type in any clinical landscape.



Last detection technologies for high technical performance



High resolution image, up to 72,5µm of voxel size



Fast acquisition



SMAR-T option, the innovative Selective Metal Artifacts Reduction Tool

Meet the Prime Family

X-MIND® prime is available in four models, 2D wall mounted, 3D wall mounted, 2D ceph and 3D ceph, to suit any installation setup and clinical need.



X-MIND® prime 2D ceph



X-MIND® prime 2D wall mounted



X-MIND® prime 3D wall mounted



X-MIND® prime 3D ceph

Save up to 30% on CBCT & OPG While stocks last

For more Information or to place orders, call **Acteon ANZ**

Toll Free **1300 ACTEON** (1300 228 366) or Sydney (02) 9669 2292 or Email acteon.anz.promo@acteongroup.com

Servotome II A New Impulse In Electrosurgery

2 YEARS WARRANTY




Experience the benefits of eletrosurgery

High frequency electrosurge or orthodontics, periodontics and other oral surgery with a single footswitch, a box with 10 electrodes & a conductive bracelet.


Electrosection

Vaporizing the tissue with the active electrode makes it possible to cut quickly without damaging the edges of the incision.





Fulguration

Makes it possible to surface coagulate gingival tissue while protecting the underlying layer.



Coagulation

Concentrating energy at the surface of a massive electrode makes it possible to diffuse the heat to the surrounding tissue and produce instantaneous hemostasis.



PIEZOELECTRIC SURGERY IS POWERFUL ENOUGH TO CARRY OUT EXTRACTIONS AT THE SAME SPEED¹ AS ROTARY INSTRUMENTS!

- Socket preservation⁹
- Safe for soft tissues
- Recommended for immediate implant placement
- No force needed
- 50% less pain & swelling²
- 98% less analgesics intake¹

1, 2, 9 References available on file



PIEZOTOME CUBE



BUY 1 x PIEZOTOME® CUBE (F50102),
GET 1 x Extraction II Tip Kit (F87546) or
 1 x Universal Tip Kit (F87539) **FREE**

Acteon Ultrasonic Cleaners

- HYGIENIC Aluminum control panel
- SILENT PEHD cover
- Stainless steel Structure / 2 years guarantee
- Waterproof switch

ULTRASONICS CLEANER 5L - 280x180x260mm
 - Usable volume 3.2L (PH500N)

ULTRASONICS CLEANER 5L WITH HEATER
 20-60°C, 280x180x260mm
 - Usable volume 3.2L with Heater (PH500CN)



Contact your preferred Acteon dealer for Special Ultrasonic Cleaner Pricing

For more Information or to place orders, call **Acteon ANZ**
 Toll Free **1300 ACTEON** (1300 228 366) or Sydney (02) 9669 2292 or Email acteon.anz.promo@acteongroup.com

SOFTLY

Amalgamator



Softly is very stable, extremely quiet, and very easy to clean and maintain. To guarantee totally safe operation, Softly is deactivated when you open the front cover to access the capsule. Softly is made entirely of injection-molded ABS material.

It complies with all electrical and mechanical safety standards.

Contact your preferred Acteon dealer for Special Softly Pricing

EXPASYL™

EXACT

An EXACT solution for atraumatic gingival retraction


Expasyl™ Exact paste provides excellent benefits with little risks. Its viscosity has been carefully calculated to open the sulcus effectively without damaging the epithelial attachment.

Your procedure is both safe and effective.

No need to wait for the tissue healing since it has not been damaged.



Generates pressure 1.7 to 9.2 times greater*



Note: Hemostasyl is included in this study but is not a retraction compound.

Adapted from American Association for Dental Research

*than other retraction paste available on the market.⁴ Data on file.

BUY 1 x Boxe 20 or 50 Pack GET 1 x Expasyl 5 Pack FREE

All clinical applications are easy with the all-in-one air polisher

AIR N GO

EASY



- A single handpiece to cover supra-gingival prophylaxis and subgingival treatments
- Greater freedom of movement for accurate treatments
- 4 interchangeable nozzles available: SUPRA, PERIO easy, PERIO, PERIO Maintenance –
- Easier cleaning of the handpiece
- The Air-N-Go Easy handpiece is suitable for SUPRA and SUB gingival treatments.
- Sodium bicarbonate-based **CLASSIC** powder
- Glycine-based **PERIO** powders



BUY 1 x AIR N GO, GET 4 x Bottles 250g AIR-N-GO Powder FREE

(Bonus bottles contain four flavours: peppermint, lemon, raspberry & cola. o/n F10255)


For more Information or to place orders, call **Acteon ANZ**

Toll Free **1300 ACTEON** (1300 228 366) or Sydney (02) 9669 2292 or Email acteon.anz.promo@acteongroup.com

RisKontrol® your 'must have' for safe, aseptic procedures

METAL TIPS

It is impossible to clean the inside of reusable tips adequately due to the presence of biofilm.⁽¹⁾



(1,2) Data on file

"With regard to cleaning and sterilization, disposable tips were found to be more acceptable in clinical use than metal tips."⁽¹⁾

Save time & money in your daily basis practice

- Single-use
- Avoid cleaning and sterilization steps

Ensure your clinical protection 20 times less risk means more safety and serenity on a daily basis, for safe care and greater patient satisfaction.

Protect your patient Single-use = reduced risk of direct or cross-contamination

- Minimise the spread of viruses



BUY 4 x Boxes Riskontrol tips GET 1 x Adaptor or 1 x Box Riskontrol FREE



NEWTRON® (Satelec) Handpieces


1. NEWTRON handpiece
Autoclavable Handpiece: GII 2-Tone (Grey) NEWTRON OEM Ref# **F12281**

2. NEWTRON® SLIM B.LED handpiece
with blue LED ring; Compatible with NEWTRON B.LED and NEWTRON XS B.LED devices for use with FLAG fluorescent liquid assistance guide Ref# **F12900**

3. NEWTRON Slim B.LED Handpiece
with white LED ring for NEWTRON B.LED, XS B.LED Ref# **F12905 3.**

NEWTRON® Tips - superior ultrasonics for all your procedures

NEWTRON® tips conserve teeth and meet all clinical requirements. NEWTRON tips were thought by dentists for dentists. They will offer you specific solutions, comfort and visibility.



Contact your Acteon Dealer for Special Pricing...

For more Information or to place orders, call **Acteon ANZ**
Toll Free **1300 ACTEON** (1300 228 366) or Sydney (02) 9669 2292 or Email acteon.anz.promo@acteongroup.com

Euronda E10 Autoclave



Digital Traceability

In compliance with EN 13060, all sterilization cycles are saved on USB (for E10) or SD Card (for E9) and can be transferred to a PC via Ethernet or WiFi connection



Load validation and release

This function allows users to keep a digital record of the compliance of the sterilization process in both the loading phase and the unloading phase. The operator ensures that the pouches are dry and checks their integrity and the state of the instruments, while also making sure that the sterilization parameters are satisfactory.

Print set options

The expandable traceability system allows three types of printers to be installed, depending on requirements. Set 1 and 2 can be installed at the front of the autoclave without any additional space whilst print set 3 is installed outside the device.



E-Memory is the software for the complete archiving, management and traceability of all sterilization cycles and data directly on a PC. This application can automatically, semi-automatically or manually capture data from the sterilization process and save it directly to the computer. This data can be easily organised, searched, printed and emailed. E-Memory also makes it possible to trace and select all information related to the cycles performed by the different devices (thermodisinfector, sealer and autoclave) involved in the sterilization process. It is easy to consult and manage, and allows you to send selected cycles and data directly to technical support in just a few clicks, and to automatically produce reports and PDF certificates with header, graphic logo and personal data.

BUY 1 x E10 Euronda Autoclave, GET 1 x in built Printer FREE
Contact your Acteon Dealer for Special Introductory Pricing...

For more Information or to place orders, contact you preferred Acteon distributor or call **Acteon ANZ**
Toll Free **1300 ACTEON** (1300 228 366) or Sydney (02) 9669 2292 or Email acteon.anz.promo@acteongroup.com

Complete the order form and email to acteon.anz.promo@acteongroup.com

- | | |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p><input type="checkbox"/> BUY 1 x Plug N Sensor
GET 1 x 1 x Sopix Plug In Dock FREE</p> <p><input type="checkbox"/> BUY 1 x Mini Led Ortho, Supercharge or Standard
GET 1 x Light Guide FREE</p> <p><input type="checkbox"/> BUY 1 x <input type="checkbox"/> 20 or <input type="checkbox"/> 50 Expasyl Exact Capsules
GET 1 x (5) Pack of Expasyl Exact FREE</p> <p><input type="checkbox"/> BUY 1 x Air-N-Go
GET 4 x Bottles 250g 4 Flavors AIR-N-GO Powder
(Peppermint, lemon, Raspberry, Cola) (Ref# F10255) FREE</p> | <p><input type="checkbox"/> BUY 1 x 20 Expasyl Capsules
GET 1 x (5) Pack of Expasyl Exact FREE</p> <p><input type="checkbox"/> BUY 4 x Boxes Riskontrol
GET 1 x <input type="checkbox"/> Box or 1 x <input type="checkbox"/> Adaptor FREE</p> <p><input type="checkbox"/> BUY 1 x Piezotome Cube (Ref# F50108)
GET 1 x <input type="checkbox"/> Extraction Tip kit II (Ref# F87546) or
<input type="checkbox"/> Universal Kit (NEW KIT) FREE</p> <p><input type="checkbox"/> BUY 1 x E8 Autoclave
GET 1 x Print Set FREE</p> |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|



New Product Release

Contact your authorised Acteon dealer for more information on Special Introductory pricing for the New Euronda E10 Autoclave.

Please print **PROMO CODE: Con N2025-1** legibly on your Invoice. To process your free good(s), email and attach a PDF of your proof of purchase to acteon.anz.promo@acteongroup.com

* The required purchase from an authorized Acteon dealer must be reflected on one invoice dated between January 1st to March 31st 2025. Invoices received after April 25th 2025 will not be honoured. Please, print promo code CON N2025-1 con legibly on your invoice. To process your free good(s), the doctor will email and attach a PDF of your proof of purchase to acteon.anz.promo@acteongroup.com Multiple redemptions on a single invoice are valid; Offer cannot be combined with any other Acteon or dealer special. Allow 4 weeks for delivery of free goods. Incomplete or illegible coupons or invoices will not be honoured. Pictures for illustration purposes only.

Name: _____

Office Name: _____

Account #: _____

Address: _____

Postcode: _____

Phone: _____

Email: _____

I would like to receive promotional, technical, and product information via email. Please send me an email so that I can opt-in to Acteon eNews.*



Acteon Australia

Suite 1, 13B Narabang Way
Belrose NSW 2085

For more Information or to place orders, contact you preferred Acteon distributor or call **Acteon ANZ**
Toll Free **1300 ACTEON** (1300 228 366) or Sydney (02) 9669 2292 or Email acteon.anz.promo@acteongroup.com