



# MATRIX

DENTAL SERVICES

To Order Tel: 1800 220 211 Fax: 9738 0822  
E-mail: sales@matrixdental.com.au



## Coricama Promotions

**Freight FREE\***

\*Orders over \$300.00



# CORICAMA

since 1873

MADE IN ITALY



OPEN



The Exclusive GAMMAFIX Locking System, comfortable in use keeps the tray perfectly closed on all conditions.

CLOSED



The mark shows the opening position (lateral) and closing position (vertical)



Open corners allow residuals and water to flow out from tray and reduce drying time



**C911030 GAMMAFIX Tray UN Red**

(185 x 73 x 35mm)

5 Instrument Position

**Regular price \$110.73**

**Promotional Price \$88.00**



C110125 Mouth Mirror Heads 4 Plane  
C110130 Mouth Mirror Heads 5 Plane  
Pack of 12

**Regular price \$28.21**

**Promotional Price \$22.80**

Silver reflection layer is at the glass-underside. The image is reflected from below the glass surface




Front Surface  
Rhodium reflection  
Layer



C110220 Mouth Mirror Heads 4 F/surface  
C110225 Mouth Mirror Heads 5 F/surface  
Pack of 12

**Regular price \$53.28**

**Promotional Price \$42.00**



**C110010 Diagnostic Kit  
1st Examination**

**Contains**

- 1 x GAMMAFIX Tray 911030
- 1 x Mouth mirror head with handle
- 1 x Explorer 23/17A
- 1 x Perio Probe CP15
- 1 x Meriam tweezer 160mm

**Regular price \$202.25**

**Promotional Price \$168.00**

## Mouth Mirror Handles



C113500  
Smooth 8mm

**Reg price \$16.90**  
**Promo Price \$12.65**



C113800  
Octagonal

**Reg price \$10.30**  
**Promo Price \$8.25**



C113860  
Hexagonal  
with scale

**Reg price \$15.50**  
**Promo Price \$11.60**



C113880  
Octagonal

**Reg price \$29.57**  
**Promo Price \$22.00**



# MATRIX

DENTAL SERVICES

A.B.N. 38 273 488 179  
Unit 4 130-132 Bayfield Road North Bayswater Vic. 3153  
Visit our Web Site matrixdental.com.au

## TWEEZERS



C161560  
Flagg  
150mm

Reg price \$17.16  
Promo Price \$12.65



C161320  
Perry  
130mm

Reg price \$17.16  
Promo Price \$12.65



C161360  
Perry  
150mm

Reg price \$17.16  
Promo Price \$12.65



C161565  
College  
150mm

Reg price \$17.16  
Promo Price \$12.65



C161770  
Meriam  
160mm

Reg price \$17.16  
Promo Price \$12.65

## EXPLORERS



C170415  
Single End  
#23  
Shepherd's hook

Reg price \$12.56  
Promo Price \$9.46




C178420  
Double End  
#23/17

Reg price \$19.80  
Promo Price \$14.85



C178430  
Double End  
#23/17A

Reg price \$19.80  
Promo Price \$14.85



C178440  
Double End  
#23/6

Reg price \$19.80  
Promo Price \$14.85



C178750  
Double End  
#6/3

Reg price \$19.80  
Promo Price \$14.85

## PROBES

C178490  
Double End  
Nabers # 2N

Reg price \$30.34  
Promo Price \$22.00



To measure horizontal  
and vertical pocket depth  
of multi-rooted teeth in  
furcation areas



C178485  
Double End  
Nabers # 2N  
Laser Markings  
3-6-9-12mm

Reg price \$37.11  
Promo Price \$27.50



C212440  
Single Ended  
CP15 UNC  
Laser Markings  
1-2-3-4-5-6-7-  
8-9-10-11-12-  
13-14-15mm

Reg price \$27.72  
Promo Price \$20.79



C212430  
Single Ended  
Laser Markings  
CP12  
3-6-9-12mm

Reg price \$27.72  
Promo Price \$20.79



C212450  
Single Ended  
WHO Ball Tip  
Laser Markings  
3.5-5.5-8.5-11.5mm

Reg price \$30.28  
Promo Price \$22.60



C214450  
Double Ended  
Williams/23  
Grooved Markings  
1-2-3-5-7-8-9-10mm

Reg price \$27.46  
Promo Price \$20.57




C214435  
Double Ended  
Goldman-Fox/23  
Grooved Markings  
1-2-3-5-7-8-9-10mm

Reg price \$27.46  
Promo Price \$20.57




C214460  
Double Ended  
Goldman-Fox  
Williams  
Round/Flat  
Grooved Markings  
1-3-5-7-8-9-10mm

Reg price \$34.06  
Promo Price \$25.55



C214485  
Double Ended  
CP12/23  
Laser Markings  
3-6-9-12mm

Reg price \$32.87  
Promo Price \$24.65



C214480  
Double Ended  
CP11/23  
Laser Markings  
3-6-8-11mm

Reg price \$32.87  
Promo Price \$24.65

Coricama Probes Single or Double Ended  
6mm Solid Handle

