



DentaPure[®] Independent Water Bottle Cartridge



Installs into your units existing water bottle pickup tube via included quick-connect fitting

- 1 installation provides 365 days (or 240L of water usage if records are kept) of safe, compliant dental unit water.
- Can use either tap or distilled water.
- Once installed, no monitoring/shocking protocol is required.
- Cartridge reduces bacterial exposure by eluting non-allergenic, elemental iodine into dental water.
- Change Every 365 days or 240L of water

Introductory price \$ 380.00

**MADE IN
USA** 