



May-June 2025 Dealer Promotion

Unit 4, 130-132 Bayfield Road East, Bayswater North, Victoria 3153 Australia
1800 220 211 | sales@matrixdental.com.au | www.matrixdental.com.au

Promotion dates: 1st May to 30th June 2025

LEDERMIX™

A Long History as the Leading Clinically Proven Solution.

Buy a Ledermix Paste and
receive the below
BONUS BUNDLE FREE!

- 1 X LEDERMIX CEMENT
- 1 X LEDERMIX HARDENER
- 1 X LEDERMIX DROPPER



**20%
OFF**

Now Only
\$159.00*
Incl. GST

**DON'T LET THEIR
PAIN BECOME
YOUR PAIN.**



**BUY 5
GET 1
FREE**

**MGuard
Mixing Tips**



Light Body Materials
Yellow #4

71004
BAG - 50

\$17⁵⁰



Regular Body Pink

71006
BAG - 50

\$17⁵⁰



C & B Material 10:1
& 4:1

71007
Bag - 50

\$14⁵⁰



Machine Mix Red for
New Penta

71013
BAG - 50

\$30⁰⁰



Intraoral Yellow
#15

71015
BAG - 100

\$13⁰⁰



Machine Mix Yellow

71014
BAG - 50

\$58⁰⁰


**Hogenspitze
Anaesthetic Needles**

BOX-100

\$13⁰⁰

- 71015 27G S 21mm Imp Thread
- 71023 27G 25mm Imp Thread
- 71021 27G L 35mm Imp Thread
- 71019 27G L 40mm Imp Thread

- 71022 30G XS 12mm Imp Thread
- 71016 30G S 21mm Imp Thread
- 71020 30G L 25mm Imp Thread



**X-ray Tabs
Self Adhesive**

Z0050
BOX-500

\$24⁰⁰

**BUY 5
GET 1
FREE**



**MGuard
Mask**



Earloop



Tie On

- Meets ASTM Standard F2100-11
- Barrier Level 2 Protection
- BFE Bacteria Filtration Efficiency 3 micron > 95%
- PFE Particle Filtration Efficiency 0.1 micron > 98%
- Listed on the TGA

- Meets ASTM Standard F2100-11
- Barrier Level 2 Protection
- Anti-Fog Filter
- BFE Bacteria Filtration Efficiency 3 microns > 95%
- PFE Particle Filtration Efficiency 0.1 micron > 98%
- Listed on TGA

- Meets ASTM Standard F2100-11
- Barrier Level 3 Protection
- Anti-Fog Filter
- BFE Bacteria Filtration Efficiency 3 microns > 95%
- BFE Particle Filtration Efficiency 0.1 micron > 98%
- Listed on TGA

Tie On Level 2
Barrier **\$6⁵⁰**
70200

Anti-Fog Tie On
Level 2 Barrier **\$7⁰⁰**
70206

Surgical Anti-Fog
Earloop Level 3 Barrier **\$13⁰⁰**
70200

Anti-Fog Earloop
Level 2 Barrier **\$6⁰⁰**
70200

**Richmond
Cotton Rolls
Braided
10mm #2**

200204
CTN - 2000

\$89⁰⁰

**BUY 3
GET 1
FREE**



**Richmond
Cotton Rolls
Economy**

BOX - 1000

#1 8mm **\$18⁰⁰**
39014

#2 10mm **\$15⁰⁰**
39015





SSWhite Bur Carbide Surgical HP

\$73⁵⁰

5 PACK Sterile

Bur Carbide Surg HP Round 018 51mm
310249 Sterile

55119



ROUND

Surgical Length (HP)	310249
Shape #	06
Head Length (mm)	1.8
Ozdent Code	55119

Bur Carbide Surg HP Round 023 51mm
310250 Sterile

55123



ROUND

Surgical Length (HP)	310250
Shape #	08
Head Length (mm)	2.3
Ozdent Code	55123

Bur Carbide Surg HP Round 023 65mm
310261 Sterile

55120



ROUND

Shank #3 (65mm) HP	30036
Shape #	08
Head Length (mm)	2.3
Ozdent Code	55120

Bur Carbide Surg HP TF 702 016 51mm
310252 Sterile

55116




TAPERED
FISSURE

Surgical Length (HP)	310252
Shape #	702
Head Length (mm)	1.6
Ozdent Code	55116

Bur Carbide Surg HP TF 702 016 65mm
310262 Sterile

55117



TAPERED
FISSURE

Shank #3 (65mm) HP	30037
Shape #	702
Head Length (mm)	1.6
Ozdent Code	55117

Bur Carbide Surg HP TF 702L 016 51mm
30024 Sterile

55125



TAPERED
FISSURE

Shank #2 (51mm) HP	30024
Shape #	702L
Head Length (mm)	1.6
Ozdent Code	55125

Bur Carbide Surg HP TF 703 021 51mm
30025 Sterile

55126



TAPERED
FISSURE

Shank #2 (51mm) HP	30025
Shape #	703
Head Length (mm)	2.1
Ozdent Code	55126

Bur Carbide Surg HP TFR 021 51mm
310256 Sterile

55115



TAPERED
FISSURE
ROUND

Shank #2 (51mm) HP	30028
Shape #	1703L
Head Length (mm)	2.1
Ozdent Code	55115

Bur Carbide Surg HP TFR 1703 021
51mm 310256 Sterile

55128



TAPERED
FISSURE
ROUND

Shank #2 (51mm) HP	30027
Shape #	1703
Head Length (mm)	2.1
Ozdent Code	55128

Bur Carbide Surg HP TFR 703L 51mm
310255 Sterile

55121



TAPERED
FISSURE
ROUND

Shank #2 (51mm) HP	30026
Shape #	703L
Head Length (mm)	2.1
Ozdent Code	55121

SSWhite Bur Carbide Surgical FG/LS

\$73⁵⁰

10 PACK Non-Sterile



Bur Carbide Surg L/S FG 1702 016
350013

55025



TAPERED
FISSURE
ROUND

Shank #2	350013
Shape #	1702
Head Length (mm)	1.6
Ozdent Code	55025

Bur Carbide Surg L/S FG TF 702 350011
55024

55024




TAPERED
FISSURE

Surgical Length (FG)	350011
Shape #	702
Head Length (mm)	1.6
Ozdent Code	55024

Bur Carbide Surg L/S FG TFR 1703
021 350014

55026



TAPERED
FISSURE

Shank #2	350014
Shape #	1703
Head Length (mm)	2.1
Ozdent Code	55026

Bur Carbide Surg L/S RA Rnd 023
55023

55023



ROUND

Surgical Length (RA)	14118-5
Shape #	08
Head Length (mm)	2.3
Ozdent Code	55023

Bur Carbide Surg L/S RA TF 702 (016)
350030

55016



TAPERED
FISSURE

Surgical Length (RA)	350030
Shape #	702
Head Length (mm)	1.6
Ozdent Code	55016

Duralast Bur Steel Round

\$36⁵⁰

25 PACK

PER PACK OF 25

- High quality slow speed steel burs
- TGA listed

61006 Bur Steel Round #1/2 006 RA

61008 Bur Steel Round #1 008 RA

61010 Bur Steel Round #2 010 RA

61012 Bur Steel Round #3 012 RA

61014 Bur Steel Round #4 014 RA

61016 Bur Steel Round #5 016 RA

61018 Bur Steel Round #6018 RA

61021 Bur Steel Round #7 021 RA

61023 Bur Steel Round #8 023 RA

61025 Bur Steel Round #9 025 RA

61027 Bur Steel Round #10 027 RA


PAYING TOO MUCH FOR STERILISATION MONITORING?

**BUY 10
GET 2
FREE**

**MGuard
Steri Pouches Self-Seal**

Internal and External
Class 4 Indicators

- Wide Triple Seals**
More secure construction less chance of instrument breaching
- Blue Tinted Laminate**
Identifies any tears or punctures in the film
- No printing within the sealing area**
No resistance to steam penetration or ink contamination
- Pre-fold**
Ensures accurate, fast seal



Corner Tasks
Prevents edge roll-up and dust contamination

Medical Grade Paper
As per ISO 11607-1

MGuard Steri Pouches Self Seal 57x100 1057100M BX - 200 \$6⁵⁰	MGuard Steri Pouches Self Seal 90x135 1090135M BX - 200 \$10⁰⁰	MGuard Steri Pouches Self Seal 70x230 1070230M BX - 200 \$11⁵⁰	MGuard Steri Pouches Self Seal 90x230 1090230M BX - 200 \$11⁰⁰	MGuard Steri Pouches Self Seal 135x260 1135260M BX - 200 \$23⁰⁰	MGuard Steri Pouches Self Seal 190x330 1190330M BX - 200 \$45⁰⁰	MGuard Steri Pouches Self Seal 250x380 1250380M BX - 200 \$45⁰⁰
--	---	---	---	--	--	--

Steam Monitor
SafeCheck Class 6 134C 3.5 min
54080 250 PACK

\$49⁰⁰

Steam Monitor
Bowie & Dick Test Pack Disposable
54090 BX - 10

\$49⁰⁰

Steam Monitor
Helix Test Kit 1x PCD + 250 strips
54085 KIT

\$189⁰⁰

**BUY 5
GET 1
FREE**

MGuard
Gown MGuard SMS Blue AAMI
Level 1 Protection
70110
10 PACK

- Infection Control gown with neck and waist ties and cotton stretch cuff
- Offers AAMI Level 1 Protection
- Cool and comfortable to wear
- TGA Listed (334546)

\$23⁰⁰

**BUY 10
GET 1
FREE**

MGuard
Face Shield w Anti-Fog
Coating
70218
EACH

- Full Length Face Shield
- Anti-Fog, Anti-Static and Anti-Glare
- Optically clear
- Thick foam head band provides generous space between face and shield
- Ideal when wearing most loupes
- Listed on the TGA (331423)

\$2⁵⁰