

# THE NEW Venus<sup>®</sup>

THE BEAUTY OF HAVING THE CHOICE

## Venus<sup>®</sup> PURE

The beauty of simplicity



## Venus<sup>®</sup> ART

The beauty of versatility



Why not choose simplicity without compromising aesthetics?

Venus responded to your demand for reduced complexity in direct restorations by dividing the brand into two categories: the new Venus Pure for simple, everyday applications and Venus Art for highly aesthetic restorations.



**KULZER**  
MITSUI CHEMICALS GROUP

# Venus<sup>®</sup> PURE

For monochromatic layering to maximise simplicity.



## Venus Pure Shades

Universal light-curing nanohybrid composite for incremental layering

Venus Pearl Pure Shades for everyday anterior and posterior restorations. The simple selection of four shades (LIGHT, MEDIUM, DARK, BLEACH) not only simplifies the procedure, but also achieves excellent results in monochromatic layering, while reducing stress and errors in shade selection.

Venus Pure Shades come in two consistencies, Venus Diamond (firm consistency) and Venus Pearl (creamier consistency).



## Venus Diamond Flow

Flowable light-curing nanohybrid composite for incremental layering

Venus Diamond Flow Pure Shades are suitable for both anterior and posterior everyday restorations. By choosing from just four shades – LIGHT, MEDIUM, DARK and BLEACH – you can streamline your process and achieve satisfactory results. In addition, this minimised shade selection reduces the likelihood of selection errors.



## Venus Bulk Flow ONE

Flowable bulk-filling composite for increments of up to 4 mm

Venus Bulk Flow ONE is the first flowable one-shade bulk-filling composite that does not require an additional capping layer. Its flowable consistency makes it the composite of choice for bulk-fillings with up to 4 mm increments in the posterior region.



# Venus<sup>®</sup> ART

For polychromatic layering to maximise aesthetics.



## Venus Art Shades

Universal light-curing nanohybrid composite for incremental layering

Venus Pearl Art Shades are for cases with exceptionally high aesthetic demands and which do not allow any compromise on craftsmanship and attention to detail. Venus Pearl Art Shades are the answer for those who want to go beyond the simple and explore the limits of what is aesthetically possible.

Venus Art Shades come in two consistencies, Venus Diamond (firm consistency) and Venus Pearl (creamier consistency).



## Venus Diamond Flow

Flowable light-curing nanohybrid composite for incremental layering

The flowable nanohybrid composite for aesthetic perfection is called Venus Diamond Flow. Among other indications, it is ideal for class V restorations. Its excellent radiopacity ensures a reliable diagnosis.



## Venus Color

Light-curing characterising material

Whether cracks in the enamel, discoloured fissures or white spots – your restorations will look entirely natural. With the five Venus Color shades, which imitate the natural appearance, you can give your patients a characteristic tooth colouration.



# Venus<sup>®</sup> PURE

## Venus Pearl Pure Shades

Shades	Syringe	PLT
	1 x 3g	20 x 0.2g
LIGHT	6609 6880	6609 6876
MEDIUM	6609 6881	6609 6877
DARK	6609 6882	6609 6878
BLEACH	6609 6883	6609 6879
ONE	6608 1836	6608 1837

## Venus Diamond Pure Shades

Shades	Syringe	PLT
	1 x 4g	20 x 0.25g
LIGHT	6609 6872	6609 6868
MEDIUM	6609 6873	6609 6869
DARK	6609 6874	6609 6870
BLEACH	6609 6875	6609 6871
ONE	6608 1838	6608 1839

## Venus Diamond Flow Pure Shades

Shades	Syringe	PLT
	1 x 1.8g	20 x 0.2g
LIGHT	6610 0268	6610 0272
MEDIUM	6610 0269	6610 0273
DARK	6610 0270	6610 0274
BLEACH	6610 0271	6610 0275

## Venus Bulk Flow ONE Refills

Content	Syringe
1x2g syringe; 5 Cannulas	6609 4723
	PLT
20x0.2g PLTs	6609 4724

## Venus Bulk Flow ONE – Value Packs

Content	Syringe
10x2g syringes; 50 Cannulas; Pictorial card	6609 5335
	PLT
80x0.2g PLTs; Pictorial card	6609 5388

## Venus Pearl Pure Kits

Content	PLT
<b>VENUS PEARL PURE KIT</b>	
20x0.2g PLTs each of Venus Pearl shades LIGHT + MEDIUM; 2g syringe of Venus Bulk Flow ONE; 2ml bottle of iBOND Universal; 5 Flow Cannulas; Translation card; Pictorial card iBOND Universal	6609 8266
<b>VENUS PEARL PURE COLLECTION KIT</b>	
20x0.2g PLTs each of Venus Pearl shades LIGHT + MEDIUM; 2g syringe of Venus Bulk Flow ONE; 1.8g syringes each of Venus Diamond Flow shades LIGHT + MEDIUM; 3x5 Flow Cannulas; Translation card	6610 0279
<b>VENUS PEARL PURE KIT</b>	Syringe
3g syringes each of Venus Pearl shades LIGHT + MEDIUM; 2g syringe of Venus Bulk Flow ONE; 2ml bottle of iBOND Universal; 5 Flow Cannulas; Translation card; Pictorial card iBOND Universal	6609 8267
<b>VENUS PEARL PURE COLLECTION KIT</b>	Syringe
3g syringes each of Venus Pearl shades LIGHT + MEDIUM; 2g syringe of Venus Bulk Flow ONE; 1.8g syringes each of Venus Diamond Flow shades LIGHT + MEDIUM; 3x5 Flow Cannulas; Translation card	6610 0278

## Venus Diamond Pure Kits

Content	PLT
<b>VENUS DIAMOND PURE KIT</b>	
20x0.25g PLTs each of Venus Diamond shades LIGHT + MEDIUM; 2g syringe of Venus Bulk Flow ONE; 2ml bottle of iBOND Universal; 5 Flow Cannulas; Translation card; Pictorial card iBOND Universal	6609 8264
<b>VENUS DIAMOND PURE COLLECTION KIT</b>	
20x0.25g PLTs each of Venus Diamond shades LIGHT + MEDIUM; 2g syringe of Venus Bulk Flow ONE; 1.8g syringes each of Venus Diamond Flow shades LIGHT + MEDIUM; 3x5 Flow Cannulas; Translation card	6610 0277
<b>VENUS DIAMOND PURE KIT</b>	Syringe
4g syringes each of Venus Diamond shades LIGHT + MEDIUM; 2g syringe of Venus Bulk Flow ONE; 2ml bottle of iBOND Universal; 5 Flow Cannulas; Translation card; Pictorial card iBOND Universal	6609 8265
<b>VENUS DIAMOND PURE COLLECTION KIT</b>	Syringe
4g syringes each of Venus Diamond shades LIGHT + MEDIUM; 2g syringe of Venus Bulk Flow ONE; 1.8g syringes each of Venus Diamond Flow shades LIGHT + MEDIUM; 3x5 Flow Cannulas; Translation card	6610 0276

## iBOND Universal

The universal, light-curing adhesive for any bonding technique and all indications. Being compatible with all dental materials, it is the all-rounder for reliable bond strength.



[kulzer.com/ibond](http://kulzer.com/ibond)

# Venus<sup>®</sup> ART

## Venus Pearl Art Shades

Shades	Syringe	PLT	PLT
	1x3g	20x0.2g	10x0.2g
<b>Universal shades</b>			
A1	6604 8175	6604 8090	
A2	6604 8176	6604 8141	
A3	6604 8178	6604 8143	
A3.5	6604 8179	6604 8144	
A4	6604 8180	6604 8145	
B1	6604 8182	6604 8147	
B2	6604 8183	6604 8148	
B3	6604 8184	6604 8149	
C1	6604 8185	6604 8150	
C2	6604 8186	6604 8161	
C3	6604 8187	6604 8162	
D3	6604 8188	6604 8163	
HKA2.5	6604 8177	6604 8142	
HKA5	6604 8181	6604 8146	
<b>Opaque dentine shades</b>			
OLC	6604 8189		6604 8164
OMC	6604 8190		6604 8165
ODC	6604 8191		6604 8093
OXDC	6604 8192		6604 8094
<b>Incisal shades</b>			
CL	6604 8196		6604 8098
AM	6604 8197		6604 8099
CO	6604 8198		6604 8100
YO	6604 8199		6604 8171
<b>Bleach and special shades</b>			
BL	6604 8193		6604 8096
BXL	6604 8194		6604 8097
OB	6604 8195		6604 8095
CORE	6604 8200		6604 8172*
GUM	6604 8201		6604 8173

\*10x0.3g

## Venus Pearl Master Assortment

Content	PLT
0.2 g x 10 PLTs (Shades A1, A2, A3, A3.5, A4, HKA2.5, ONE); 0.2 g x 5 PLTs (Shades B1, B2, B3, C1, C2, C3, D3, CL, AM, CO, YO, OLC, OMC, ODC, OXDC, BL, BXL, OB, GUM, HKA5); 0.3 g x 5 PLTs (Shade CORE); Venus Pearl Layering Guide	6609 6867

## Venus Color Art Shades

Content	Syringe	Syringe
	1x1g	5x1g
<b>Shades</b>		
WHITE	6605 4206	
BLUE	6605 4207	
CORN	6605 4208	
AMBER	6605 4209	
CHOCO	6605 4210	
<b>Assortment</b>		
Venus Color Syringe Assortment (1x1g each: white, blue, corn, amber, choco)		6605 4221

## Hygienic Protective Sleeve

Content	
Hygienic protective sleeves to protect syringes from contamination by blood, saliva or other contaminants. 300 pieces	6609 7644

## Venus Diamond Art Shades

Shades	Syringe	PLT	PLT
	1x4g	20x0.25g	10x0.25g
<b>Universal shades</b>			
A1	6603 5471	6603 9008	
A2	6603 5472	6603 9009	
A3	6603 5474	6603 9011	
A3.5	6603 5475	6603 9012	
A4	6603 5476		6603 9013
B1	6603 5478	6603 9015	
B2	6603 5479	6603 9016	
B3	6603 5480		6603 9017
C2	6603 5481	6603 9018	
C3	6603 5482		6603 9019
D3	6603 5483		6603 9020
HKA2.5	6603 5473	6603 9010	
HKA5	6603 5477		6603 9014
<b>Opaque dentine shades</b>			
OL	6603 5484		6603 9021
OLC	6604 6735		6604 6731
OM	6603 5485	6603 9022	
OMC	6604 6736		6604 6732
OD	6603 5486	6603 9023	
ODC	6604 6737		6604 6733
OXDC	6604 6738		6604 6734
<b>Incisal shades</b>			
CL	6603 5490		6603 9027
AM	6603 5491		6603 9028
CO	6603 5492		6603 9029
YO	6603 5493		6603 9030
<b>Bleach shades</b>			
BL	6603 5487		6603 9025
BXL	6603 5488		6603 9026
OB	6603 5489		6603 9024

## Venus Diamond Flow Art Shades

Shades	Syringe	PLT
	1x1.8g	20x0.2g
<b>Universal shades</b>		
A1	6604 0354	6604 0368
A2	6604 0355	6604 0369
A3	6604 0358	6604 0371
A3.5	6604 0359	6604 0372
A4	6604 0360	6604 0373
B2	6604 0362	6604 0375
HKA2.5	6604 0357	6604 0370
HKA5	6604 0361	6604 0374
<b>Opaque dentine shades</b>		
OM	6604 0363	6604 0376
Baseliner	6604 1052	6604 1051
<b>Incisal shades</b>		
CL	6604 0364	6604 0377
<b>Bleach shades</b>		
BXL	6604 0365	6604 0378

## Shade Guides

Content	
Venus Diamond Shade Guide	6603 9004
Venus Pearl Shade Guide	6605 5463




# UNLEASHING STRENGTH: THE UNIQUE ADVANTAGES OF OUR TCD MATRIX

The TCD matrix features the TCD urethane molecule, which vibrates less and takes up less space than conventional molecules like Bis-GMA/TEGDMA. This limited movement brings notable advantages, such as lower shrinkage stress, reducing the likelihood of marginal gaps and secondary cavities. Additionally, the extensive cross-linking in the TCD matrix enhances the composite's strength and longevity, thereby lowering the risk of chipping or breaking.

Moreover, the TCD matrix is entirely free of BPA structures. Testing by an external accredited laboratory has confirmed that both unpolymerised and polymerised Venus Pearl and Venus Diamond composites show no detectable BPA release.<sup>1</sup>

Mean shrinkage stress (blue squares) after 24 h [MPa] and mean flexural strength (bars) after 24 h [MPa] of the tested composite materials



<sup>1</sup>Source: Test report 10-2023

<sup>2</sup>Ceram-x Spectra® ST HV is a trademark of Dentsply Sirona Inc. · Filtek® Universal is a trademark of 3M Company · Tetric® Prime is a trademark of Ivoclar Vivadent AG · Omnichroma® is a trademark of Tokuyama Corporation · Essentia® is a registered trademark of GC Corporation · Admira® Fusion is a registered trademark of Voco GmbH

Source: Flexural strength and polymerisation shrinkage stress of composite materials, Leyer A, Schweppe J, Utterodt A, J Dent Res Vol 103 (Spec Iss B): 0026, 2024

# BE PART OF IT!

JOIN OUR COMMUNITY

SHARE YOUR CASES & DISCOVER  
AN ACTIVE COMMUNITY THAT LOVES  
SMILES AS MUCH AS WE DO.



kulzer.com/  
venus-bepartofit