



Disposable Dental Needles

Ultra sharp needles
for patients comfort

Triple-bevel lancet tip

Bevel tip indicator for
ease of guidance.

Flexible cannula for
reduced risk of breakage

Sterile **STERILE EO**

Latex Free



- Triple-bevel lancet point without coring for easy tissue penetration



- Hot sealed plastic hub & paper sealed label



Code	Brand	Description	Qty/Ctn
D-DN2738	DIAGURU	27 G x 38 mm - Green 100Pcs/box	50 boxes
D-DN3021	DIAGURU	30 G x 21 mm - Yellow 100Pcs/box	50 boxes



Contact our Customer Service
to place your order

Matrix Dental
U4/130-132 Bayfield Rd. East,
Bayswater North, VIC, AUSTRALIA

Phone 1800 220 211
sales@matrixdental.com.au