



# HOGEN spitze

## Dental Needles

Hogenspitz Dental Cartridge Needles with tri-bevel point and silicon coated for smooth tissue penetration, bevel indicating mark on the hub, designed to be used with imperial and metric threaded syringes.



Dental Needles With Bevel Indicator  
On The Hub **27G Short 21mm**

Code - Z1015  
Box - 100



Dental Needles With Bevel Indicator  
On The Hub **30G Short 21mm**

Code - Z1016  
Box - 100



Dental Needles With Bevel Indicator  
On The Hub **27G XX Long 40mm**

Code - Z1019  
Box - 100



Dental Needles With Bevel Indicator  
On The Hub **30G Long 25mm**

Code - Z1020  
Box - 100

