

Kerr™

Exclusive Offers Through



Discount
Up to
36%
Off

Valid While Stocks Last

Elements IC Obturation Dual System

~~\$6,281.55~~

\$4,838.90



DA DENTAL ADVISOR

Winner of the 2023 Top Award for
ENDODONTIC OBTURATION SYSTEM!

MXKERTV01

Elements IC Downpack System

~~\$3,312.93~~

\$2,310.00



MXKERTV02

Say
"GOODBYE"
to your Bunsen Burner

All prices are inclusive of GST

Powered obturation with no strings attached.

Obturation is a critical step in the root canal procedure as it helps to prevent leakage, entomb remaining bacteria and reduce the chance of reinfection. The elements IC obturation system was designed to be used with the Continuous Wave warm vertical condensation technique. The new extended battery life and Inductive Charging Base makes the elements IC the latest technology for 3D endodontic obturation without cords or wires getting in the way.



Features & Benefits | A System Based Approach.



Cordless emulation of the Elements Obturation Unit
 » same Continuous Wave™ downpack tips
 » same backfilling cartridges
 » same performance-without the cord



Ergonomic design for comfortable and balanced grip



Batteries fully charge in less than 3 hours



Inductive Charging:
 NO ELECTRICAL CONNECTOR
 between handpieces and base




Gutta-Percha cartridge level indicator,
 for better planning of the procedure



360° Ring on both handpieces for control and comfort


Downpack Device



- DOUBLE battery run time compared to elements free, more than 16 canals on a single charge
- Priming time at least 20% shorter compared to elements free
- Compatible with most common cleaning solutions (CaviWipes, Sani-Cloth Plus, Optim 33TB)
- Adjustable temperature range of 100°C - 230°C

Indicator	4	3	2	1	0
Gutta Percha Level	Full	75%	50%	25%	0%

The indicator estimates how many canals you can obturate with the material left in the GP cartridge.



Backfill Device

Ordering Information

Part #	Item Description
973-0600-TYPEI	elements IC Obturation System
973-0602-TYPEI	Downpack Unit
973-0604-TYPEI	Backfill Unit
973-0610	Dual Charger
973-0612	Single Charger
973-0615	Transformer
973-0616-TYPEI	Power Cord
973-0620	Battery



Heat Pluggers

Part #Item Description

952-0004	Fine Buchanan Plugger, .06 Taper
952-0005	Fine-Medium Buchanan Plugger, .08 Taper
952-0006	Medium Buchanan Plugger, .10 Taper
952-0007	Medium-Large Buchanan Plugger, .12 Taper
952-0031	Extra-fine Buchanan Plugger, .04 Taper



Hand Pluggers

Part # Item Description

974-0058	Buchanan Hand Plugger-1-Red
974-0059	Buchanan Hand Plugger-2-Blue
974-0060	Buchanan Hand Plugger-0-Yellow



elements Cartridge

Part #	Item Description
972-1002	Gutta Percha, Medium Body, 23 Gauge - Silver
972-1005	Gutta Percha, Heavy Body, 23 Gauge - Silver
972-1003	Gutta Percha, Light Body, 25 Gauge - Silver
972-2500	Gutta Percha, Medium Body, 23 Gauge - Gold
972-2502	Gutta Percha, Heavy Body, 23 Gauge - Gold
972-2501	Gutta Percha, Light Body, 25 Gauge - Gold

