



AeroAcoustic blade

Reverse flow protection

Refined ceramic bearings



Glass rod fiber optics



Automatic pressure regulation 1



Scratch resistance Iconic Prime body



Powerful & Precise Power Head - 25W



Incredibly Quiet Only 55 dB noise



Small Head Coupling **Power Head** Connection 25 Watt 18 Watt HP21KL HP22KL Kavo NSK HP21NL HP22NL W&H HP21WL HP22WL Sirona HP21SL - N/A -

\$1,683.00



Save \$500 with Trade-In handpiece of

- ✓ Any brand
- ✓ Any type
- ✓ Any condition



Now with 2 Years FairCare Warranty²

Unconditional warranty except physical damage on Prime Line High-speeds

V COUPLINGS





 QC6016KW \$495.00

Non-optic QC4014K

\$296.00



 QC6016NT \$385.00

Non-optic

 QC4014NT \$296.00



with fibre optic

 QC5016W \$385.00

Non-optic

 QC4014W \$296.00



Sirona connection with LED & water

 QC6016SW \$495.00

Non-optic -NA-

1) Only available for Kavo coupling 2) FairCare warranty applicable on Prime Line high speeds sold after 1st Oct 2018

- **Buy 3 Get 15% OFF**
- Buy 4 Get 20% OFF

Mix & Match





2 Years FairCare Warranty **Unconditional** warranty except physical damage









- % 1:5 Red •
- O Max. 200,000 rpm
- O Internal water
- For FG-burs (Ø 1.6 mm)
- LP15L (F/Optic) RRP: \$1,826.00

Buy 3 - \$1,552.10 ea

Buy 4 - \$1,460.80 ea

- % 1:1 Blue •
- C Max. 40,000 rpm
- O Internal water
- For CA-burs (Ø 2.35 mm)
- LP11L (F/Optic) RRP: \$1,287.00

Buy 3 - \$1,093.95 ea

Buy 4 - \$1,029.60 ea

- % 5.4:1 Green •
- C Max. 7,407 rpm
- O Internal water
- For CA-burs (Ø 2.35 mm)
- LP54L (F/Optic) RRP: \$1,320.00

Buy 3 - \$1,122.00 ea

Buy4 - \$1,056.00 ea

- % 1:1 Blue ●
- C Max. 40,000 rpm
- O Internal water
- For HP-burs (Ø 2.35 mm)
- LP01L (F/Optic) RRP: \$1,045.00

Buy 3 - \$888.25 ea

Buy 4 - \$836.00 ea

AIR MOTORS W





- Internal water
- LED
- AM1006 \$990.00



- ☑ Internal water
- AM1004 \$748.00







AM0004 \$473.00

NEW

KAVO COMPATIBLE HEADS V



- Kavo 68LH compatible
- Blue (1:1)
- LCH11 \$ 440.00



- Kavo 67LH compatible
- Green (2:1)
- LCH21 \$ 440.00



- Kavo 31I R compatible
- Green (2:1)
- LCH21P \$ 440.00





C



62 dB Operating Noise



AeroPerformance blade



Refined ceramic bearings



Glass rod fiber optics



Smooth Eco body



• \$890.00 (F/Optic)

Available for KaVo*, NSK*, W&H* and Sirona* Couplings



BUY 3 or MORE Only

3 x Eco Line HS with F/O for any coupling connection

> Savings from \$392

Coupling	Power Head	Small Head	Power Head	Small Head
Connection	24W (F/Optic)	16W (F/Optic)	24W (N/Optic)	16W (N/Optic)
KaVo	HE21KL	HE22KL	HE21K	HE22K
NSK	HE21NL	HE22NL	HE21N	HE22N
W&H	HE21WL	HE22WL	HE21W	HE22W
Sirona	HE21SL	– N/A –	– N/A –	– N/A –
• \$ 890 00 (E/Optic)			• \$649.00 (N/Optic)	

ECO LINE POWER TURBINE





28 Watt Power



Buv 3 - Get 10% OFF

60 dB Operating Noise



AeroPerformance blade



Glass rod fiber optics



Refined ceramic bearings



Smooth Eco body

HE20KL \$1,265.00



Available for only KaVo Multiflex coupling fitting

ECO LINE FIXED CONNECTION TURBINES

- Direct 4 hole connection
- Refine ceramic bearings
- Push Button System
- HE17 Standard head with 4 hole water spray
- HE18 Small head with 3 hole water spray



Buy 3 - Get 15% OFF Mix & Match



- % 1:5 Red O Internal water
- For FG-burs (Ø 1.6 mm)
- LE15L (F/Optic) RRP: \$1,485.00 Buy 3 - \$1,262.25 e
- LE15 (Non-optic) RRP: \$1,287.00 Buy 3 - \$1,093.95 ea

% 1:1 Blue •

- O Internal water
- For CA-burs (Ø 2.35 mm)
- LE11L (F/Optic) RRP: \$990.00
- Buy 3 \$ 841.50 e
- LE11 (Non-optic) RRP:\$770.00

Buy 3 - \$ 654.50 ea



FairCare Warranty

year unconditional warranty except physical damage



- % 5.4:1 Green •
- O Internal water
- For Screw-in and Snap-on Prophy Cups
- LE54LP (F/Optic) RRP: \$1,210.00
- Buy 3 \$1,028.50 ea
- LE54P (Non-optic) RRP: \$935.00

Buy 3 - \$ 794.75 ea



- % 1:1 Blue •
- O Internal water
- For HP-burs (Ø 2.35 mm)
- LE01L (F/Optic) RRP: \$935.00
- Buy 3 \$ 794.75 ea
- LE01 (Non-optic) RRP: \$748.00
 - Buy 3 \$ 635.80 ea

BASIC LINE CONTRA - ANGLES V



- % 1:1 Blue •
- O Push button • LB11/11
- Buy 1 \$418.00
- Buy 3 \$349.00 ea



- % 1:1 Blue O Latch Type
- LB11/11C
- Buy 1 \$418.00

Buy 3 - \$349.00 ea



- 4:1 Green•
- O Snap-on Cups
- LB41/11P
- Buy 1 \$ 462.00

Buy 3 - \$396.00 ea



- % 1:1 Blue •
- LB02
- Buy 1 \$ 418.00

Buy 3 - \$349.00 ea



Buv 3 - Get 15% OF





MK-deut PROPHY LINE



EMS Non-optic Type EN-041

 SC21E Buy 1 - \$462.00

Buy 3 - \$392.70 ea

EMS Piezon Non-optic Type EN-061

 SC21E2 Buy 1 - \$561.00 Buv 3 - \$476.85 ea



EMS Piezon LED Type EN-060A

 SC21EL Buy 1 - \$748.00 Buv 3 - \$635.80 ea



KaVo Piezo LED Type 1.007.3995

 SC21KL* NEW Buy 1 - \$748.00

Buy 3 - \$635.80 ea *Only for use with EMS compatible scaler tips

Type F 12 281 SC21S Buy 1 - \$ 462.00

Acteon-Satelec Non-optic

Buy 3 - \$392.70 ea



Acteon-Satelec LED Type F 12 609

 SC21SL Buy 1 - \$748.00

Buy 3 - \$635.80 ea

V SCALER TIPS



- High-quality Diamond-Like-Carbon coating
- Extremely hard and wear-resistant surface

• SCALING DLC COATED Only available for Satelec fitting

ST02SD SATELEC 1S \$88.00

ST03SD **EMS A** \$88.00

ST04SD SATALEC 2 \$132.00

PERIO DLC COATED Only available for Satelec fitting

ST06SD SATELEC 10Z \$132.00

ST10SD SATALEC 10P \$154.00

ST08SD SATALEC S07 \$154.00

ENDO DLC

SCALING STAINLESS STEEL Only available for Satelec fitting

ST02S SATELEC 1S \$63.80

ST03S **EMSA** \$63.80

ST04S SATELEC 2 \$88.00

PERIO STAINLESS STEEL Available for EMS & Satelec fittings

ST07E/ST07S **EMS DPL3** SATELEC 10X \$110.00

ST09E/ST09S **EMS HPL**

\$110.00 \$110.00

ST10E/ST10S SATALEC 10P

ST11E/ST11S EMS RT1 \$110.00

EMS and Acteon-Satelec are registered trademarks. MK-dent GmbH has no economic connection with the aforementioned companies.

Universal Torque Wrench

TOZQ CONTZOL®

Precision instrument with a Torque Adjustment System for tightening of Prosthetic components at specific Torque, for any implant system.



Automatic Declutching when Selected Torque is reached

- Compliance with Protocols for screw tightening
- Limits the risk of the screw fracturing or loosening



Universal connection of **Prosthetic mandrels**

- 7 different Torque levels (from 10 to 35 N.cm)
- Can be used with **ALL Implant Systems**



Micro-Head and Ergonomic Handle

 Easy access in restricted area with Excellent Visibility



Easy and Quick to use

 Reduced motion compared with using a ratchet wrench while tightening




Lightweight ONLY 135 gms

Order code: 15501

3 Easy Steps









www.anthogyr.com




V PROPHY LINE AIR POLISHING DEVICE



Buy 1 Polisher and receive 1 bottle of Prophy Powder FREE and 50% discount on 1 additional tip

• PR1011 \$1,375.00



air-driven polishing device

 360° rotatable outer handle due to unique RotaGrip® for simple and fatique-proof operation



- Prophy Powder (300 gms) \$60.50 each
- R Lemon PP1012

Additional Prophy tip



Coupling connection order codes:

Kavo - PR1011K, NSK - PR1011N, W&H - PR1011W, Sirona - PR1011S

• **SP3011** \$374.00

V LEDs



Buy 4 - Get 10% OFF

- Quality illumination, No Glare
- 5500K Natural Daylight Colour
- Clear Vision & Sharper Contrast
- Longer Durability & Economical



\$88.00

BU8012 - KaVo & MK-dent coupling

BU8012KM - KaVo & MK-dent motors

BU8012WH - W&H handpieces
BU8012WLT - W&H airmotors LT25

BU8012WEA - W&H motors 40LT, 50 LT & A51LT

BU8012SG - Sirona coupling and motors

BU8012N - NSK coupling

BU8012NT - NSK motors TiM40 & Ti205L

BU8012NM - NSK coupling KCL,SCL,PTL LED

♥ PREMIUM LUBRICATION OIL



- Fully Synthetic
- For any brand of handpieces and motors
- LU1011 (500 ml) \$68.20



- Fully Synthetic
- For KaVo
 Quattrocare 2

 lubrication system
- LU1022 (500 ml) \$ 68.20



- For W&H Assistina lubrication system
- LU1031 (500 ml) \$88.00

MANUFACTURER



MK-dent GmbH Marie-Curie-Str. 2 22941 Bargteheide GERMANY

www.mk-dent.com

YOUR SPECIALIST DEALER