

# FEBRUARY - APRIL 2024 SPECIAL!

TOP OF THE CLASS IN PAEDIATRIC CROWNS

## ACERO CROWNS<sup>3S</sup>

**ACERO CROWNS<sup>3S</sup>** crowns are crafted from surgical grade 316L stainless steel, which is more resistant to corrosion than the 304 stainless steel used by competitors. Our crowns are preformed and pre-contoured (commonly referred to as "pre-crimped") with a shallow occlusal anatomy to preserve healthy tooth structure with less tooth reduction. The bell-shaped crowns are manufactured with optimal thickness, and they replicate the contacts and embrasures of the natural dentition. The laser-etched crowns with easy-to-read numbers won't rub off after sterilisation but can be removed with polishing points if desired. **ACERO CROWNS<sup>3S</sup>** crowns are proudly made in the USA and are 100% interchangeable with other leading brands but at a competitive price.

Crown S/S Acero 3S

Primary Kit-48

Code: 22001

Kit-48

\$449.<sup>00</sup>

Crown S/S Acero 3S

Primary Kit-240

Code: 22000

Kit-240

\$2150.<sup>00</sup>

INTRO  
PRICE

### Advantages:

- Competitive cost against other leading brands.
- Laser-etched sizes (autoclave proof = less waste).
- Surgical grade 316L US stainless steel.
- 100% compatible and interchangeable with the leading brands.
- Preformed and pre-contoured to prevent micro-gaps.
- Made in the USA.



INTRODUCING

## ACERO NARROW<sup>NA</sup>

**ACERO NARROW** crowns were developed to address the challenges of restoring primary first molars with large interproximal caries. **ACERO NARROW** crowns maintain all the desirable traits of our popular molar crowns but with a 12% reduction in mesial width to compensate for the space loss due to mesial drift of the 2nd primary molar.

Crown S/S Acero Narrow

Primary Molar Kit-120


Code: 22200

Kit-120

\$999.<sup>00</sup>

INTRO  
PRICE

REGULAR ACERO NARROW



### Advantages:

- Less mesial distal width to treat that large class 2 caries space lost.
- Acero case features lightweight design with easy closure and no spillage out of the compartments.



## THE MOST DURABLE PAEDIATRIC CROWN ON THE MARKET

### ACERO CROWNS<sup>ES</sup>



**ACERO CROWNS<sup>ES</sup>** Aesthetic Anterior Crowns are manufactured using 316L surgical-grade stainless steel veneered with durable white resin and are proven to be highly fracture-resistant compared to other brands. **ACERO CROWNS<sup>ES</sup>** are available in a variety of sizes and created with optimal thickness to replicate the natural dentition. Our crowns are labelled using the **card on inner roof of the case depicts Australian system** with letters A through T to identify the corresponding primary tooth with ease.

Crown S/S Resin Coated

Acero ES Tray Kit-30

Code: 22130

Kit-30

**\$1649.00**

**INTRO  
PRICE**

#### Advantages:

- Less tooth destruction required to place ES Crowns.
- Preformed and pre-contoured.
- Life-like translucency.
- Acero case features lightweight design with easy closure and no spillage out of the compartments.



INTRODUCING

### ACERO CROWNS<sup>Z1</sup>



Continuing our mission to bring the highest quality paediatric crowns as well as giving dentists the most value, **ACERO** is proud to launch our new **Z**, preformed zirconia crowns. The **Z Crowns** are simply beautiful with the morphology and translucency of the natural tooth. Our molar crowns have a flat occlusal inner surface with retentive boxes making tooth preparation simple and quick while providing great retention for the crowns.

Crown Zirconium Acero Z1

Anteriors Tray Kit-48

Code: 22250

Kit-48

**\$2849.00**

**INTRO  
PRICE**

#### Advantages:

- Features strong, biocompatible and aesthetic qualities. Prized in dentistry for years.
- Acero case features lightweight design with easy closure and no spillage out of the compartments.



# FEBRUARY - APRIL 2024 SPECIAL!

Bur Carbide Fissurotomy  
Original 310215  
Code: 55140  
Pkt-5



**10% OFF!**

\$41.<sup>50</sup>

Bur Carbide Fissurotomy  
NTF 310218  
Code: 55141  
Pkt-5



**10% OFF!**

\$41.<sup>50</sup>

Bur Carbide Fissurotomy  
STF 18013  
Code: 55142  
Pkt-5



**10% OFF!**

\$41.<sup>50</sup>

Forcep Extraction Duralast  
Milk Tooth Tweezer  
Code: 27497  
Each



**20% OFF!**

\$67.<sup>00</sup>

Forcep Extraction Duralast  
#37 Child Upper Incisor  
Code: 27437  
Each



**20% OFF!**

\$110.<sup>00</sup>

Forcep Extraction Duralast  
#37N Child Upper Anterior  
Code: 27438  
Each



**20% OFF!**

\$110.<sup>00</sup>

Forcep Extraction Duralast  
#39 Child Upper Molar  
Code: 27489  
Each



**20% OFF!**

\$110.<sup>00</sup>

Forcep Extraction Duralast  
#123 Child Lower/Root Uni  
Code: 27492  
Each



**20% OFF!**

\$110.<sup>00</sup>

Forcep Extraction Duralast  
#138 Child Upper Teeth  
Code: 27493  
Each



**20% OFF!**

\$110.<sup>00</sup>

Forcep Extraction Duralast  
#158 Child Upper Teeth  
Code: 27499  
Each



**20% OFF!**

\$110.<sup>00</sup>

Forcep Extraction Duralast  
#159 Child Upper Teeth  
Code: 27496  
Each



**20% OFF!**

\$110.<sup>00</sup>

Forcep Extraction Duralast  
#160 Child Lower  
Molar Hawks  
Code: 27490  
Each



**20% OFF!**

\$110.<sup>00</sup>

Forcep Extraction  
Duralast 161 Child  
Lower Molar Hawk  
Code: 27491  
Each



**20% OFF!**

\$110.<sup>00</sup>

Forcep Extraction  
Duralast #162 Child  
Lower Anterior  
Code: 27495  
Each



**20% OFF!**

\$110.<sup>00</sup>

Mouth Prop Child  
Pedo Silicon  
Non-Latex Yellow  
Code: 60210  
Pkt-2



**BUY 2  
GET 1 PKT-2  
ADULT  
MOUTH PROP  
BONUS!**

\$15.<sup>50</sup>

Mouth Prop Child  
Silicon Non-Latex  
Green  
Code: 60212  
Pkt-2



**BUY 2  
GET 1 PKT-2  
ADULT  
MOUTH PROP  
BONUS!**

\$15.<sup>50</sup>

TePe Select Compact  
Soft Toothbrush  
Code: 84872  
Each



**30% OFF!**

\$3.<sup>50</sup>

TePe Select Compact  
X-Soft 'Zoo'  
Illustrations Cel  
Code: 84610  
Each



**30% OFF!**

\$4.<sup>50</sup>

TePe Select Compact  
X-Soft Toothbrush  
Blisterpack  
Code: 84838  
Each



**30% OFF!**

\$4.<sup>50</sup>

TePe Select Compact  
X-Soft Toothbrush  
Code: 84873  
Each



**30% OFF!**

\$3.<sup>50</sup>

TePe Select Medium  
Toothbrush Blisterpack  
Code: 84831  
Each



**30% OFF!**

\$3.<sup>50</sup>

TePe Select Mini  
X-Soft Toothbrush  
Code: 84163  
Each



**30% OFF!**

\$2.<sup>50</sup>

Brush Bur Stainless  
Steel Bristle  
Code: 46100  
Each



**20% OFF!**

\$4.<sup>00</sup>

Brush Instrument  
Cleaning Nylon Bristle  
Code: 46105  
Each



**20% OFF!**

\$4.<sup>00</sup>

Varnish Alpha Pro  
NAF Melon Export  
Code: 20070  
Box-50



**20% OFF!**

\$115.<sup>00</sup>

UVEX Pheos Child  
White/Orange Arm  
HC-AF Tinted Lens  
Code: 9192201  
Each



**BUY 2  
GET 1 UVEX  
ADULT SPEC  
BONUS!**

\$29.<sup>50</sup>

UVEX Pheos Child  
White/Green Arm  
HC-AF Clear Lens  
Code: 9192200  
Each



**BUY 2  
GET 1 UVEX  
ADULT SPEC  
BONUS!**

\$29.<sup>50</sup>

UVEX Pheos Small  
Foam Guard  
Code: 9192002  
Each



**BUY 2  
GET 1 UVEX  
ADULT SPEC  
BONUS!**

\$7.<sup>50</sup>



## FEBRUARY - APRIL 2024 SPECIAL!

Sterilisation Reel  
(5cm x 200m)  
Code: 70300  
Each

**20% OFF!**

**\$89.00**



Sterilisation Reel  
(7.5cm x 200m)  
Code: 70302  
Each

**20% OFF!**

**\$109.00**



Sterilisation Reel  
(10cm x 200m)  
Code: 70304  
Each

**20% OFF!**

**\$139.00**



Sterilisation Reel  
(20cm x 200m)  
Code: 70307  
Each

**20% OFF!**

**\$259.00**



Steam Monitor SafeCheck  
Class 6 134°C (3.5 min)  
Code: 54080  
Pkt-250

**BUY 3  
GET 1  
BONUS!**


**\$55.00**



Steam Monitor Helix Test  
Kit 1 x PCD + 250 Strips  
Code: 54085  
Kit

**BUY 2  
GET 1  
BONUS!**

**\$259.00**



Steam Monitor Bowie &  
Dick Test Pack Disposable  
Code: 54090  
Kit-20

**20% OFF!**

**\$99.00**




### Self Seal Pouches With Class 4 Indicators

#### INDICATIONS

**Wide Triple Seals Plus  
Extra Chevron Seal**  
More secure construction  
less chance of  
instrument breaching.

**Green Tinted Laminate**  
Identifies any tears or  
punctures in the film.

**No printing within  
the sealing area**  
No resistance to  
steam penetration  
or ink contamination.



**Corner Tasks**  
Prevent edge  
roll-up & dust  
contamination.

**Two Internal  
Class 1 Indicators**  
Release strip +  
preprinted provides  
positive proof  
that steam has  
penetrated the  
pouch.

**Medical  
Grade Paper**  
As per ISO  
11607-1.

**Pre-fold**  
Ensure accurate,  
fast seal.

**BUY 3  
GET 1  
BONUS!**



#### FEATURES

- Compliance With ISO Standard 1
- Strong seals 10mm width
- Safe print design
- Pack traceability
- Medical grade paper as per Standard 11607 - 1
- Sufficient glue overlap
- Internal indicators
- Pre-folded

Steri Pouches  
Self Seal MGuard  
(57mm x 100mm)  
Code: 1057100M  
Box-200

**\$6.50**

Steri Pouches  
Self Seal MGuard  
(90mm x 135mm)  
Code: 1090135M  
Box-200

**\$9.00**

Steri Pouches  
Self Seal MGuard  
(70mm x 230mm)  
Code: 1070230M  
Box-200

**\$10.00**

Steri Pouches  
Self Seal MGuard  
(90mm x 230mm)  
Code: 1090230M  
Box-200

**\$12.00**

Steri Pouches  
Self Seal MGuard  
(135mm x 260mm)  
Code: 1135260M  
Box-200

**\$19.00**

Steri Pouches  
Self Seal MGuard  
(190mm x 330mm)  
Code: 1190330M  
Box-200

**\$33.00**

Steri Pouches  
Self Seal MGuard  
(250mm x 380mm)  
Code: 1250380M  
Box-200

**\$45.00**